

## Quick Shop Plus

1st Edition  
Spring 2011

### 100-Well Freezer Boxes

- Polypropylene, autoclavable
- Store down to -90°C
- Black well indexes
- Colored, hinged lid




ORDER #	DESCRIPTION	QTY	SALE!
FRB-100	100-Well Freezer Box, Natural	5/pk	\$20.00
FRB-100C	100-Well Freezer Box, Assorted*	5/pk	\$24.00
FRB-100(x)	100-Well Freezer Box, Specify Color**	5/pk	\$23.00

\* Blue, Green, Pink, Yellow, Orange

\*\* Replace (x) with B (Blue), G (Green), P (Pink), Y (Yellow), O (Orange)

### 81-Well Freezer Boxes

- Polypropylene, autoclavable
- Store down to -90°C
- Black well indexes
- Clear, lift-off lid



ORDER #	DESCRIPTION	QTY	SALE!
FB-81	81-Well Freezer Box, Natural	5/pk	\$20.00
FB-81C	81-Well Freezer Box, Assorted*	5/pk	\$24.00
FB-81(x)	81-Well Freezer Box, Specify Color**	5/pk	\$23.00

\* Blue, Green, Pink, Yellow, Orange

\*\* Replace (x) with B (Blue), G (Green), P (Pink), Y (Yellow), O (Orange)

### Polycarbonate Cryostorage Boxes

- Store down to -196°C
- Numerical index stamped on box base
- Drain holes in boxes and lids



ORDER #	DESCRIPTION	QTY	SALE!
CF-25	25-Well PC Cryostorage Box, Assorted*	5/pk	\$19.00
CF-25(x)	25-Well PC Cryostorage Box, Specify Color**	5/pk	\$19.00
CF-81	81-Well PC Cryostorage Box, Assorted*	5/pk	\$29.00
CF-81(x)	81-Well PC Cryostorage Box, Specify Color**	5/pk	\$29.00
CF-100	100-Well PC Cryostorage Box, Assorted*	5/pk	\$31.00
CF-100(x)	100-Well PC Cryostorage Box, Specify Color**	5/pk	\$31.00

\* Blue, Green, Pink, Yellow, Orange

\*\* Replace (x) with B (Blue), G (Green), R (Red), Y (Yellow), O (Orange)


### SlipLock Flipper® Interlocking Racks

- Use individually, interlocked, or stacked
- 24 x 1.5/2ml and 16 x 0.5ml on one side
- 60 x 0.2ml on the other side
- Supplied as 4 racks with 2 lids



ORDER #	DESCRIPTION	QTY	SALE!
SLF-RACK	SlipLock Flipper Rack, 4 Racks w/ 2 Lids	1/pk	\$27.00


## Revolver Rack Workstation



- 4 separate, interlocking, different color racks
- Holds multiple 0.2ml tubes and PCR plates, 0.5ml, 1.5/2.0ml, 15ml and 50ml tubes

ORDER #	DESCRIPTION	QTY	SALE!
REV-RACK	Revolver Rack, Assorted, 4-piece	1/pk	\$30.00

## 96-Well PCR Rack



- 8 x 12 format, 0.2ml
- Holds tubes, strip tubes, and plates (except full skirt plates)

ORDER #	DESCRIPTION	QTY	SALE!
WR-96A	96-Well PCR Rack, Assorted*	5/pk	\$23.00
WR-96(x)	96-Well PCR Rack, Specify Color**	5/pk	\$23.00

\* Blue, Green, Pink, Yellow, Orange  
\*\* Replace (x) with B (Blue), G (Green), P (Pink), Y (Yellow), O (Orange)

## 96-Well PCR Work-Up Rack



- 8 x 12 format, 0.2ml
- Stackable
- Holds tubes, strip tubes, and plates (including full skirt plates)

ORDER #	DESCRIPTION	QTY	SALE!
WUR-96A	96-Well PCR Work-Up Rack, Assorted*	5/pk	\$23.00

\* Blue, Green, Pink, Yellow, Orange

## 96-Well Flipper® Rack



- Holds 96 1.5/2.0ml or 0.5ml Tubes
- Polypropylene, autoclavable

ORDER #	DESCRIPTION	QTY	SALE!
FL-96-05	96 Flipper Rack w/ Lid, Natural	5/pk	\$31.00
FL-96F-5R	96 Flipper Rack w/ Lid, Assorted*	5/pk	\$32.00
FL-96F(x)-5R	96 Flipper Rack w/ Lid, Specify Color**	5/pk	\$32.00

\* Blue, Green, Pink, Yellow, Orange  
\*\* Replace (x) with B (Blue), G (Green), P (Pink), Y (Yellow), O (Orange)

## 80-Well Microtube Racks



- Polypropylene, autoclavable
- Frosted writing area

ORDER #	DESCRIPTION	QTY	SALE!
36010	80-Well Microtube Rack, Natural	5/pk	\$15.00
3601F	80-Well Microtube Rack, Assorted*	5/pk	\$16.00
3601(x)	80-Well Microtube Rack, Specify Color**	5/pk	\$16.00

\* Blue, Green, Pink, Yellow, Orange  
\*\* Replace (x) with B (Blue), G (Green), P (Pink), Y (Yellow), O (Orange)

## 4-Way Flipper® Racks



- Polypropylene, autoclavable
- Store down to -90°C
- Can be interlocked to expand capacity
- Holds (4) 50ml, (12) 15ml, (32) 1.5/2ml, and (32) 0.5ml tubes

ORDER #	DESCRIPTION	QTY	SALE!
FL-4N	4-Way Flipper Rack, Natural	5/pk	\$21.00
FL-4A	4-Way Flipper Rack, Assorted*	5/pk	\$23.00
FL-4(x)	4-Way Flipper Rack, Specify Color**	5/pk	\$22.00

\* Blue, Green, Pink, Yellow, Orange  
\*\* Replace (x) with B (Blue), G (Green), P (Pink), Y (Yellow), O (Orange)


## IsoFreeze<sup>®</sup> Microtube Chiller Racks



- Blue, maintains between -10°C and -20°C for up to 3 hours
- White, maintains 0°C for minimum of 5 hours
- Holds 20 x 1.5/2.0ml or 20 x 0.5ml tubes

ORDER #	DESCRIPTION	QTY	SALE!
TFL-20	Blue IsoFreeze Microtube Chiller Rack, -20°C	1/pk	\$33.00
TFL-0	White IsoFreeze Microtube Chiller Rack, 0°C	1/pk	\$53.00


## IsoFreeze<sup>®</sup> MCT Chiller Racks



- Maintains temperature less than 4°C for up to 4 hours
- Changes color with warming
- Holds 24 x 0.5/1.5/2.0ml tubes
- Supplied with lids

ORDER #	DESCRIPTION	QTY	SALE!
5630-T4	IsoFreeze MCT Rack, Purple to Pink	1/pk	\$40.00
5630-T6	IsoFreeze MCT Rack, Green to Yellow	1/pk	\$40.00

## IsoFreeze<sup>®</sup> PCR Chiller Racks



- Maintains temperature less than 4°C for up to 4 hours
- Changes color with warming
- 96-well, holds 0.2ml tubes, strips, or plates
- Supplied with lids

ORDER #	DESCRIPTION	QTY	SALE!
5640-T4	IsoFreeze PCR Racks, Purple to Pink	2/pk	\$47.00
5640-T6	IsoFreeze PCR Racks, Green to Yellow	2/pk	\$47.00

## IsoFreeze<sup>®</sup> PCR SBS Chiller Racks



- Maintains temperature less than 4°C for up to 3 hours
- Changes color with warming
- 96-well, holds 0.2ml tubes, strips, or plates
- SBS footprint for automated systems

ORDER #	DESCRIPTION	QTY	SALE!
5650-T4	IsoFreeze PCR SBS Racks, Purple to Pink	2/pk	\$49.00
5650-T6	IsoFreeze PCR SBS Racks, Green to Yellow	2/pk	\$49.00

## Floater Rack



- Can be used in water baths, ice buckets, or self-contained
- Holds 24 x 1.5/2ml and 14 x 0.5ml tubes
- Supplied as rack and hinged container

ORDER #	DESCRIPTION	QTY	SALE!
FR24	Floater Rack	1/pk	\$7.00

## 15ml and 50ml Conical Tube Freezer Boxes




- Numerically indexed
- Green, 36 x 15ml
- Blue, 16 x 50ml

ORDER #	DESCRIPTION	QTY	SALE!
5570-02	15ml Conical Tube Rack, Green	2/pk	\$20.00
5580-01	50ml Conical Tube Rack, Blue	2/pk	\$20.00

## TPP Conical Tube Rack


- Collapsible to save space
- Autoclavable
- Holds 18 x 15ml or 10 x 50ml tubes
- Alphanumerically marked for easy reference



ORDER #	DESCRIPTION	QTY	SALE!
TP80617	TPP Conical Tube Rack	1/pk	\$6.00

## CLASSIC BLUE™ Film BX

- Double emulsion
- Exceptional clarity with extremely low backgrounds
- High signal sensitivity
- Excellent lot-to-lot consistency
- Publishable grade autoradiographs
- Packaged in boxes of 100



ORDER #	DESCRIPTION	QTY	PRICE
BX57	5" x 7" Classic Blue™ Film BX	100	\$62.00
BX810	8" x 10" Classic Blue™ Film BX	100	\$84.00
BX1417	14" x 17" Classic Blue™ Film BX	100	\$199.00

## INCUSHAKE™ 10L Incubated Shaker

- Ambient +5°C to 60°C
- 30-300rpm
- Flask capacity: 5x2L, 9x1L (see website for full table)
- Supplied with non-slip rubber mat platform



ORDER #	DESCRIPTION	PRICE
IS-1010	INCUSHAKE 10L Incubated Shaker	\$4895.00
IS-1010-PLAT-F	MAGic Clamp Magnetic Platform, 18" x 18"	FREE!

Please visit our website for clamps and accessories.

## ORBI-SHAKE™ CO<sub>2</sub> Resistant Shaker

Safe for use in CO<sub>2</sub> and other extreme environments!

- Designed for use in CO<sub>2</sub> incubators without corroding
- Remote controllable – adjust settings without opening incubator door
- Speed Range: 30 to 300rpm
- Operating Temp: 0-60°C
- Change flask clamps almost instantly with MAGic Clamp magnetic platforms
- Supplied with 13" x 12" rubber-mat platform



Perfect for use with TPP Bioreactor Tubes and Bottles

Search keyword: BIOREACTOR



ORDER #	DESCRIPTION	PRICE
OS-4000	ORBI-SHAKE CO <sub>2</sub> Resistant Shaker	\$2995.00
OS-PLAT-F	MAGic Clamp Magnetic Platform, 14" x 12"	FREE!

Please visit our website for clamps and accessories.